



Product Datasheet

RPA1 Antibody (orb1563328)


Description	RPA1 Antibody
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Purified recombinant human RPA70 protein fragments expressed in E.coli.
Target	RPA1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	4 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	Calculated MW: 68 kDa; Observed MW: 68 kDa
Uniprot ID	P27694
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Expiration Date	12 months from date of receipt.



Immunoprecipitation analysis of RPA70 (8...



Western blot analysis of RPA70 in HeLa, ...



Immunocytochemistry analysis of RPA70 (8...