



Product Datasheet

NUP98 Antibody (orb1563366)

Description	NUP98 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Purified recombinant human NUP98 protein fragments expressed in E.coli.
Target	NUP98
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 198 kDa; Observed MW: 98 kDa
Uniprot ID	P52948
Dilution Range	WB: 1/500-1/1000 IP: 1/20
Expiration Date	12 months from date of receipt.



Western blot analysis of NUP98 in HeLa, ...



Immunoprecipitation analysis of NUP98 (3...