



Product Datasheet

FGFR1 Antibody (orb1563446)


Description	FGFR1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human FGFR1 around the phosphorylation site of Tyr654.
Target	FGFR1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Calculated MW: 92 kDa; Observed MW: 120 kDa
Uniprot ID	P11362
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200 ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of Phospho-FGFR1 (...)



Western blot analysis of Phospho-FGFR1 (...)



EnzymeLinked Immunosorbent Assay (Phospho...)