



Product Datasheet

Phospho-Tyrosine (6G6) Mouse mAb Antibody (orb1563492)

Description	Phospho-Tyrosine (6G6) Mouse mAb Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Purified Protein
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.



Western blot analysis of Phospho-Tyrosin...



Immunohistochemistry analysis of paraffi...