



## Product Datasheet

### DYNLT1 Antibody (orb1563611)


|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | DYNLT1 Antibody                                                                                                                         |
| <b>Species/Host</b>        | Rabbit                                                                                                                                  |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Tested Applications</b> | ICC, IHC-Fr, IHC-P, IP, WB                                                                                                              |
| <b>Immunogen</b>           | Recombinant protein of human TCTEX-1                                                                                                    |
| <b>Target</b>              | DYNLT1                                                                                                                                  |
| <b>Preservatives</b>       | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin                                              |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                  |
| <b>Concentration</b>       | 0.3 mg/ml                                                                                                                               |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>Isotype</b>             | IgG                                                                                                                                     |
| <b>Clonality</b>           | Monoclonal                                                                                                                              |
| <b>MW</b>                  | Calculated MW: 12 kDa; Observed MW: 12 kDa                                                                                              |
| <b>Uniprot ID</b>          | <b>P63172</b>                                                                                                                           |
| <b>Dilution Range</b>      | WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200<br>IP: 1/20                                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                         |



Immunocytochemistry analysis of Dynein L...



Western blot analysis of TCTEX1 in 3T3, ...



Western blot analysis of TCTEX1 in HeLa, ...