



## Product Datasheet

### PUS1 Antibody (orb1564365)

<b>Description</b>	PUS1 Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Hamster, Human, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IP, WB
<b>Immunogen</b>	Recombinant protein of human PUS1
<b>Target</b>	PUS1
<b>Preservatives</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.3 mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>MW</b>	Calculated MW: 47 kDa; Observed MW: 44,47 kDa
<b>Uniprot ID</b>	<a href="#">Q9Y606</a>
<b>Dilution Range</b>	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of PUS1 in Jurkat,...



Immunocytochemistry analysis of PUS1 (gr...