


Product Datasheet

AHSG Antibody (orb1564911)

Description	AHSG Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic peptide of human Fetuin A
Target	AHSG
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Form/Appearance	Liquid
Concentration	0.3 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	Calculated MW: 39 kDa; Observed MW: 45-50 kDa
Uniprot ID	P02765
Dilution Range	WB: 1/500-1/1000 IP: 1/20
Expiration Date	12 months from date of receipt.



Western blot analysis of Fetuin A in Hep...