



Product Datasheet

ZHX2 Antibody (orb1565453)

Description	ZHX2 Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Recombinant protein of human ZHX2
Target	ZHX2
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Form/Appearance	Liquid
Concentration	0.3 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	Calculated MW: 92 kDa; Observed MW: 110 kDa
Uniprot ID	Q9Y6X8
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200
Expiration Date	12 months from date of receipt.



Western blot analysis of ZHX2 in rat Bra...



Immunocytochemistry analysis of ZHX2 (gr...