



Product Datasheet

RHOA/RHOB/RHOC Antibody (orb1565565)


Description	RHOA/RHOB/RHOC Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr, IHC-P, WB
Immunogen	Recombinant protein of human Rho C
Target	RHOA/RHOB/RHOC
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Uniprot ID	P62745 , P61586 , P08134
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100
Expiration Date	12 months from date of receipt.



Western blot analysis of RhoA/B/C in Jur...



Immunocytochemistry analysis of Rho A/B/...



Immunohistochemistry analysis of paraffi...