



Product Datasheet

Cnot7 Antibody (orb1565570)


Description	Cnot7 Antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, IP, WB
Immunogen	A synthetic peptide of mouse CNOT7
Target	Cnot7
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Uniprot ID	Q60809
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20 FC: 1/50-1/100
Expiration Date	12 months from date of receipt.



Western blot analysis of CNOT7 in F9, C6...



Western blot analysis of CNOT7 in Molt4,...



Immunohistochemistry analysis of paraffi...