

Product Datasheet

DDX3(Phospho-Thr322) antibody (orb156560)

Description	DDX3(Phospho-Thr322) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human DDX3 around the phosphorylation site of Thr322 (VA(p-T)PG)
Target	phospho-DDX3 (Thr322)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	73 kDa
Uniprot ID	O00571
Expiration Date	12 months from date of receipt.