


Product Datasheet

TUBB Antibody (HRP) (orb1566572)

Description	TUBB Antibody (HRP)
Species/Host	Mouse
Reactivity	Hamster, Human, Monkey, Mouse, Rat
Conjugation	HRP
Tested Applications	WB
Immunogen	Purified recombinant human Tubulin beta protein fragments expressed in E.coli.
Target	TUBB
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	0.31 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 50 kDa; Observed MW: 55 kDa
Uniprot ID	P07437
Dilution Range	WB: 1/5000-1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of Tubulin beta in...