



Product Datasheet

ACTB Antibody (HRP) (orb1566573)


| | |
|----------------------------|---|
| Description | ACTB Antibody (HRP) |
| Species/Host | Mouse |
| Reactivity | Algae, Drosophila, Goat, Hamster, Human, Monkey, Mouse, Plant, Rabbit, Rat, Saccharomyces, Yeast, Zebrafish |
| Conjugation | HRP |
| Tested Applications | WB |
| Immunogen | Purified recombinant human beta-actin protein fragments expressed in E.coli. |
| Target | ACTB |
| Preservatives | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3. |
| Form/Appearance | Liquid |
| Concentration | 2.4 mg/ml |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG2b |
| Clonality | Monoclonal |
| MW | Calculated MW: 42 kDa; Observed MW: 42 kDa |
| Uniprot ID | P60709 |
| Dilution Range | WB: 1/5000-1/10000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of beta actin in H...



Western blot analysis of beta actin in r...



Western blot analysis of beta actin in r...