



Product Datasheet

TUBB Antibody (orb1566584)


Description	TUBB Antibody
Species/Host	Mouse
Reactivity	Algae, Goat, Hamster, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Purified recombinant human Tubulin beta protein fragments expressed in E.coli.
Target	TUBB
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	5.06 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 50 kDa; Observed MW: 55 kDa
Uniprot ID	P07437
Dilution Range	WB: 1/5000-1/10000 IF: 1/100-1/200
Expiration Date	12 months from date of receipt.



Western blot analysis of Tubulin beta in...



Immunocytochemistry analysis of Tubulin ...



Western blot analysis of beta Tubulin (1...