



Product Datasheet

TUBB Antibody (orb1566597)


Description	TUBB Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Synthetic Peptide of β -tubulin
Target	TUBB
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 50 kDa; Observed MW: 55 kDa
Uniprot ID	P07437
Dilution Range	WB: 1/5000-1/10000 IF: 1/100-1/200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa, mouse Bra...



Immunofluorescence analysis of beta Tubu...



Western blot analysis of β Tubulin in Ju...