

Product Datasheet

Trenbolone (orb1566843)

| | |
|----------------------------|--|
| Description | Trenbolone |
| Tested Applications | ELISA, WB |
| Preservatives | 0.2M PBS, pH7.2, 10% glycerol |
| Concentration | 1mg/ml |
| Storage | This product is stable for up to twelve months from date of receipt at 4°C, Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. |
| Note | For research use only |
| Formula | $C_{18}H_{22}O_2$ |
| Source | This compound was modified in house and tested by animal assays and ELISA. |
| CAS Number | 10161-33-8 |
| Expiration Date | 12 months from date of receipt. |