

Product Datasheet

25-OH-VitaMin D3/Calcifediol Monohydrate (orb1566847)

Description	25-OH-VitaMin D3/Calcifediol Monohydrate
Tested Applications	ELISA, WB
Preservatives	0.2M PBS, pH7.2, 10% glycerol
Concentration	1mg/ml
Storage	This product is stable for up to twelve months from date of receipt at 4°C, Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	$C_{27}H_{44}O_2$
Source	This compound was modified in house and tested by animal assays and ELISA.
CAS Number	63283-36-3
Expiration Date	12 months from date of receipt.