

## Product Datasheet

### Oxfendazole (orb1566850)

<b>Description</b>	Oxfendazole
<b>Reactivity</b>	Other
<b>Tested Applications</b>	ELISA, WB
<b>Preservatives</b>	0.2M PBS, pH7.2, 10% glycerol
<b>Concentration</b>	1mg/ml
<b>Storage</b>	This product is stable for up to twelve months from date of receipt at 4°C, Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{15}H_{13}N_3O_3S$
<b>Source</b>	This compound was modified in house and tested by animal assays and ELISA.
<b>CAS Number</b>	53716-50-0
<b>Expiration Date</b>	12 months from date of receipt.