

Product Datasheet

Streptomycin (orb1566883)

Description	Streptomycin
Reactivity	Bacteria
Tested Applications	ELISA, WB
Preservatives	0.2M PBS, pH7.2, 10% glycerol
Concentration	1mg/ml
Storage	This product is stable for up to twelve months from date of receipt at 4°C, Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	$C_{21}H_{39}N_7O_{12}$
Source	This compound was modified in house and tested by animal assays and ELISA.
CAS Number	57-92-1
Expiration Date	12 months from date of receipt.