

Product Datasheet

Ciprofloxacin (orb1566923)

Description	Ciprofloxacin
Tested Applications	ELISA, WB
Preservatives	0.2M PBS, pH7.2, 10% glycerol
Concentration	1mg/ml
Storage	This product is stable for up to twelve months from date of receipt at 4°C. Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	$C_{17}H_{18}FN_3O_3$
Source	This compound was modified in house and tested by animal assays and ELISA.
CAS Number	85721-33-1
Expiration Date	12 months from date of receipt.