

Product Datasheet

eNOS (Phospho-Ser615) antibody (orb156734)

Description	eNOS (Phospho-Ser615) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human eNOS around the phosphorylation site of Ser615 (FN(p-S)IS)
Target	Phospho-eNOS (Ser615)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	133 kDa
Uniprot ID	P29474
Expiration Date	12 months from date of receipt.