


Product Datasheet

ENTPD2 antibody (orb156742)

Description	ENTPD2 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human ENTPD2/CD39L1 (401-495/495 aa)
Target	ENTPD2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54 kDa
Uniprot ID	Q9Y5L3
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse pancreas ...