
Product Datasheet

Rabbit Anti-Sheep Cell Hemolysin (orb1568208)

Description

Rabbit Anti-Sheep Cell Hemolysin

Species/Host

Rabbit

Conjugation

Unconjugated

Tested Applications

HAI

Immunogen

Sheep washed pooled Red Blood Cells (RBC)

Preservatives

Preservative: None. Stabilizer: None. None

Form/Appearance

Lyophilized

Concentration

75.0 mg/mL

Storage

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Note

For research use only

Application notes

Rabbit Anti-Sheep Cell Hemolysin has been tested by agglutination and is suitable for agglutination, complement fixation related assays. Reconstitution Buffer: Restore with deionized water (or equivalent). Reconstitution Volume: 2.0 mL

Clonality

Polyclonal

Purity

Rabbit Sheep Cell Hemolysin was prepared from polyspecific antiserum by a delipidation and defibrination.

Expiration Date

12 months from date of receipt.