

Product Datasheet

Rabbit Anti-Sheep Cell Hemolysin (orb1568208)

Description	Rabbit Anti-Sheep Cell Hemolysin
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	HAI
Immunogen	Sheep washed pooled Red Blood Cells (RBC)
Preservatives	Preservative: None. Stabilizer: None. None
Form/Appearance	Lyophilized
Concentration	75.0 mg/mL
Storage	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	Rabbit Anti-Sheep Cell Hemolysin has been tested by agglutination and is suitable for agglutination, complement fixation related assays. Reconstitution Buffer: Restore with deionized water (or equivalent). Reconstitution Volume: 2.0 mL
Clonality	Polyclonal
Purity	Rabbit Sheep Cell Hemolysin was prepared from polyspecific antiserum by a delipidation and defibrination.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)