

Product Datasheet

Sheep Red Blood Cell RBC Antibody (orb1568259)

Description	Sheep Red Blood Cell RBC Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	HAI
Immunogen	Sheep washed pooled Red Blood Cells (RBC)
Preservatives	0.01% (w/v) Sodium Azide
Form/Appearance	Lyophilized
Concentration	85 mg/mL
Storage	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only
Application notes	Tested for agglutination of cells on titer plates. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Isotype	Antiserum
Clonality	Polyclonal
Purity	This product was prepared from polyspecific antiserum by delipidation and defibrination.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)