

Product Datasheet

Dopamine D2 receptor agonist-2 (orb1568302)

Description	Dopamine D2 receptor agonist-2 (compound 36) is a potent dopamine D2 receptor biased agonism ligand ...
Note	For research use only
Formula	$C_{25}H_{31}Cl_2N_5OS$
MW	520.52
CAS Number	[1610591-93-9]
SMILES	<chem>O=C(C1=C(N)C2=C(C)C=C(C)N=C2S1)NCCCCCN3CCN(C4=CC=CC(Cl)=C4Cl)CC3</chem>
Expiration Date	12 months from date of receipt.