

## Product Datasheet

### Azido-PEG5-acid (orb1568347)

<b>Description</b>	Azido-PEG5-acid is a non-cleavable 5 unit PEG ADC linker used in the synthesis of antibody-drug conj...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{13}H_{25}N_3O_7$
<b>MW</b>	335.3535
<b>CAS Number</b>	[1425973-16-5]
<b>SMILES</b>	<chem>O(C([H])([H])C([H])([H])OC([H])([H])C([H])([H])OC([H])([H])C([H])([H])N=[N+]=[N-])C([H])([H])C([H])([H])OC([H])([H])C([H])([H])OC([H])([H])C([H])([H])C(=O)O[H]</chem>
<b>Expiration Date</b>	12 months from date of receipt.