



Product Datasheet

CD74 Antibody (PE) (orb1570904)


Description	CD74 (the MHC II-associated invariant chain, Ii) is a type II transmembrane protein expressed in ant...
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	SU-DHL-4 lymphoma cells
Target	CD74
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 20 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P04233
Entrez	972
Expiration Date	12 months from date of receipt.



Flow cytometry multicolor surface staining...



Separation of human CD74 positive B cell...



Flow cytometry surface staining pattern ...