



Product Datasheet

CD172a Antibody (FITC) (orb1570907)

Description	CD172a, the signal-regulatory protein alpha (SIRP alpha), also known as SH2 domain-containing phosph...
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	Kg-1a cell line
Target	CD172a
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Uniprot ID	P78324
Entrez	140885
Expiration Date	12 months from date of receipt.



Separation of human monocytes (red-fille...



Flow cytometry surface staining pattern ...