
Product Datasheet

**Canine CD3 epsilon Protein, His Tag (MALS
verified) (orb1570922)**

Description

Canine CD3 epsilon Protein, His Tag

(HEK293). It contains AA Gln 2...

Reactivity

Canine

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

CD3 epsilon

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 16.7 kDa. The protein migrates as 19-25 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

Gln 22 - Leu 122

Purity

90%

Source

Canine CD3 epsilon Protein, His Tag (MALS verified) (orb1570922) is expressed from human 293 cells (HEK293). It contains AA Gln 22 - Leu 122 (Accession # P27597).

MW

16.7 kDa

Uniprot ID

P27597

Expiration Date

6 months from date of receipt.