
Product Datasheet

**Rabbit Angiopoietin-2 / ANGPT2 Protein, His Tag
(MALS verified) (orb1570923)**

Description

Rabbit Angiotensinogen Converting Enzyme 2 Protein,

human 293 cells (HEK293). It contains AA Lys...

Reactivity

Rabbit

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Angiotensin-2

Preservatives

20 mM MOPS, 5mM CHAPS, 150 mM NaCl, pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 27.4 kDa. The protein migrates as 30-34 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

Lys 275 - Phe 496

Purity

95%

Source

Rabbit Angiotensinogen Converting Enzyme 2 Protein, His Tag (orb1570923) is expressed from human 293 cells (HEK293). It contains AA Lys 275 - Phe 496 (Accession # G1T9B3-1).

MW

27.4 kDa

Uniprot ID**G1T9B3****Expiration Date**

6 months from date of receipt.