

---

## Product Datasheet

### Rabbit Angiopoietin-2 / ANGPT2 Protein, His Tag (orb1570928)

**Description**

Rabbit Angiopoietin-2 Protein, His

contains AA Tyr 19 - Phe 496 ...

**Reactivity**

Rabbit

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

Angiopoietin-2

**Preservatives**

20 mM MOPS, 150 mM NaCl, 5 mM CHAPS, pH7.5

**Form/Appearance**

Powder

**Storage**

-20°C

**Tag**

C-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 56.7 kDa. The protein migrates as 60-68 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Protein Sequence**

Tyr 19 - Phe 496

**Purity**

90%

**Source**

Rabbit Angiopoietin-2 Protein, His Tag (orb1570928) is expressed from CHO cells. It contains AA Tyr 19 - Phe 496 (Accession # G1T9B3-1).

**MW**

56.7 kDa

**Uniprot ID**
[G1T9B3](#)
**Expiration Date**

6 months from date of receipt.