
Product Datasheet

**Human LYPD1 / PHTS Protein, Mouse IgG2a Fc Tag
(MALS verified) (orb1570931)**

Description

Human LYPD1 Protein, Mouse

human 293 cells (HEK293). It contains AA L...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

LYPD1

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-mFc(mIgG2a)

Note

For research use only

Application notes

This protein carries a mouse IgG2a Fc tag at the C-terminus. The protein has a calculated MW of 37.3 kDa. The protein migrates as 43-48 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

Leu 21 - Gly 115

Purity

95%

Source

Human LYPD1 Protein, Mouse IgG2a Fc Tag (orb1570931) is expressed from human 293 cells (HEK293). It contains AA Leu 21 - Gly 115 (Accession # Q8N2G4-1).

MW

37.3 kDa

Expiration Date

6 months from date of receipt.