



## Product Datasheet

### GFPT1 antibody (orb157116)


|                            |  |
|----------------------------|--|
| <b>Description</b>         | GFPT1 antibody   |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | ELISA, ICC, IF, IHC-Fr, IHC-P, WB  |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human GFPT1 (601-699/699 aa)       |
| <b>Target</b>              | GFPT1  |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.            |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Concentration</b>       | 1mg/ml   |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>MW</b>                  | 79 kDa   |
| <b>Uniprot ID</b>          | <b>Q06210</b>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |



Immunohistochemical staining of human lu...



Western blot analysis of hela lysate usi...



Western blot analysis of 293T cell lysat...