

Product Datasheet

SARS-CoV-2 N protein Antibody (Biotin) (orb1572099)

| | |
|----------------------------|---|
| Description | SARS-CoV-2 N protein Antibody (Biotin) |
| Species/Host | Rabbit |
| Conjugation | Biotin |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant SARS-CoV-2 Nucleocapsid Protein |
| Target | SARS-CoV-2 N protein |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | ≥1mg/ml |
| Storage | Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P0DTC2 |
| Dilution Range | WB=1:500-2000,ELISA=1:5000-10000 |
| Expiration Date | 12 months from date of receipt. |