

Product Datasheet

Phospho-HSD17B13 (Ser33) Antibody (Biotin) (orb1572135)

Description	Phospho-HSD17B13 (Ser33) Antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ICC, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human HSD17B13 around the phosphorylation site of Ser33
Target	Phospho-HSD17B13 (Ser33)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7Z5P4
Dilution Range	WB=1:500-2000, ICC=1:100
Expiration Date	12 months from date of receipt.