
Product Datasheet

RNA polymerase II CTD repeat YSPTSPS (phospho S5) Antibody (BF594) (orb1572584)

Description

RNA polymerase II CTD repeat

(BF594)

Species/Host

Rabbit

Reactivity

Human

Conjugation

BF594

Tested Applications

FC, IF

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human RNA polymerase II CTD around the phosphorylation site of Ser5

Target

RNA polymerase II CTD repeat
YSPTSPS (phospho S5)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID**P24928****Dilution Range**

Flow-Cyt=1ug/Test,IF=1:100-500

Expiration Date

12 months from date of receipt.