

---

## Product Datasheet

### Phospho-Histone H3 (Thr3) Antibody (BF680) (orb1573389)

**Description**

Phospho-Histone H3 (Thr3) Antibody

|                            |   |
|----------------------------|---|
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Rat  |
| <b>Conjugation</b>         | BF680   |
| <b>Tested Applications</b> | FC, ICC, IF   |
| <b>Immunogen</b>           | KLH conjugated Synthesised phosphopeptide derived from human Histone H3 around the phosphorylation site of Thr3 |
| <b>Target</b>              | Phospho-Histone H3 (Thr3)   |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.   |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Concentration</b>       | 1mg/ml  |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.                                |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>CAS Number</b>          | Histone H3.1  |
| <b>Uniprot ID</b>          | <a href="#">P68431</a>  |
| <b>Dilution Range</b>      | Flow-Cyt=1µg /Test,IF=1:100-500,ICC=1:100-500   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |