
Product Datasheet

**Phospho-Gluccorticoid Receptor (Ser226)
Antibody (BF680) (orb1574054)**

Description

Phospho-Glucocorticoid Receptor

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF680
Tested Applications	ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Glucocorticoid Receptor around the phosphorylation site of Ser226
Target	Phospho-Glucocorticoid Receptor (Ser226)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04150
Dilution Range	IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.