

## Product Datasheet

### CD11a Antibody (BF680) (orb1574987)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | CD11a Antibody (BF680)   |
| <b>Species/Host</b>        | Mouse  |
| <b>Reactivity</b>          | Human  |
| <b>Conjugation</b>         | BF680  |
| <b>Tested Applications</b> | FC, ICC, IF  |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human CD11a                        |
| <b>Target</b>              | CD11a  |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.            |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Concentration</b>       | 1mg/ml   |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Monoclonal   |
| <b>Uniprot ID</b>          | <a href="#">P20701</a>   |
| <b>Dilution Range</b>      | Flow-Cyt=1ug/Test,IF=1:50,ICC=1:50   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |