


Product Datasheet

IRAK1 (Phospho-Thr209) antibody (orb157703)

Description	IRAK1 (Phospho-Thr209) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human IRAK1 around the phosphorylation site of Thr209 (QG(p-T)CN)
Target	Phospho-IRAK1 (Thr209)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	78 kDa
Uniprot ID	Q62406
Dilution Range	IHC-P=1:100-500, ELISA=1:5000-10000, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...