

## Product Datasheet

### KCNJ1(Phospho-Ser44) antibody (orb157734)

<b>Description</b>	KCNJ1(Phospho-Ser44) antibody
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, ICC, IF, IHC-Fr, IHC-P
<b>Immunogen</b>	KLH conjugated synthesised phosphopeptide derived from human KCNJ1 around the phosphorylation site of Ser44 (LV(p-S)KD)
<b>Target</b>	phospho-KCNJ1 (Ser44)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	45 kDa
<b>Uniprot ID</b>	<a href="#">P48048</a>
<b>Expiration Date</b>	12 months from date of receipt.