

## Product Datasheet

### Kir6.2(Phospho-Thr224) antibody (orb157754)

<b>Description</b>	Kir6.2(Phospho-Thr224) antibody
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	KLH conjugated synthesised phosphopeptide derived from human Kir6.2 around the phosphorylation site of Thr224 (KT(p-T)SP)
<b>Target</b>	phospho-Kir6.2 (Thr224)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	44 kDa
<b>Uniprot ID</b>	<a href="#">Q14654</a>
<b>Expiration Date</b>	12 months from date of receipt.