


Product Datasheet

LAT (Phospho-Tyr161) antibody (orb157776)

Description	LAT (Phospho-Tyr161) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human LAT isoform a around the phosphorylation site of Tyr161 (PG(p-Y)LV)
Target	Phospho-LAT (Tyr161)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29 kDa
Uniprot ID	043561
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...