
Product Datasheet

**Histone H3 (tri methyl K36) Antibody (BF680)
(orb1577877)**

Description

Histone H3 (tri methyl K36) Antibody

| | |
|----------------------------|--|
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | BF680 |
| Tested Applications | FC, ICC, IF |
| Immunogen | KLH conjugated synthesised methylpeptide derived from human Histone H3 around the methylation site of tri methyl K36 |
| Target | Histone H3 (tri methyl K36) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| CAS Number | Histone H3.1 |
| Uniprot ID | P68431 |
| Dilution Range | Flow-Cyt=1µg /Test,IF=1:100-500,ICC=1:100-500 |
| Expiration Date | 12 months from date of receipt. |