

Product Datasheet

Phospho-PDGFR α (Tyr754) Antibody (BF680) (orb1578191)

Description	Phospho-PDGFR α (Tyr754) Antibody (BF680)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF680
Tested Applications	ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PDGF Receptor alpha around the phosphorylation site of Tyr754
Target	Phospho-PDGFR α (Tyr754)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16234
Dilution Range	IF=1:100-500, ICC=1:100
Expiration Date	12 months from date of receipt.