

## Product Datasheet

# Phospho-PDGFR $\alpha$ (Tyr754) Antibody (BF647) (orb1578192)

<b>Description</b>	Phospho-PDGFR $\alpha$ (Tyr754) Antibody (BF647)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	BF647
<b>Tested Applications</b>	ICC, IF
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from human PDGF Receptor alpha around the phosphorylation site of Tyr754
<b>Target</b>	Phospho-PDGFR $\alpha$ (Tyr754)
<b>Preservatives</b>	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>P16234</b>
<b>Dilution Range</b>	IF=1:100-500,ICC=1:100
<b>Expiration Date</b>	12 months from date of receipt.