

Product Datasheet

Phospho-MSK1 (Ser376) Antibody (PE) (orb1578251)

Description	Phospho-MSK1 (Ser376) Antibody (PE)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	PE
Tested Applications	ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MSK1 around the phosphorylation site of Ser376 GY(p-S)FV
Target	RPS6KA5
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	90 kDa
Uniprot ID	O75582
Dilution Range	WB=1:500-2000, IHC-P=1:20-200, IHC-F=1:20-200, ICC/IF=1:100, IF=1:100
Expiration Date	12 months from date of receipt.