

Product Datasheet

MEPCE Immunizing peptide (orb15783)

| | |
|--------------------------|--|
| Description | MEPCE Immunizing peptide |
| Preservatives | 100µg of Lyophilised peptide, formulation De-ionized (DI) water. |
| Storage | Shipped at ambient temperature, store at - 20°C |
| Note | For research use only |
| Application notes | This peptide was used as immunogen to produce orb15783 (MEPCE antibody) |
| Protein Sequence | KQQRKFQYGN YCK |
| Expiration Date | 12 months from date of receipt. |