

Product Datasheet

Hsc70 Rabbit Polyclonal Antibody (PerCP) (orb1579023)

Catalog Number	orb1579023
Category	Antibodies
Description	Hsc70 Rabbit Polyclonal Antibody (PerCP)
Target	HSPA8
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	PerCP
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Bovine, Gallus, Porcine, Rabbit
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin and 0.02% Proclin300.
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated synthetic peptide derived from human HSP71 (501-600/646aa)
Tested applications	FC, ICC, IF
Dilution range	Flow-Cyt=1ug/test, ICC/IF=1:100-500, IF=1:100-500
Antibody Type	Primary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage

Store at 2-8°C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)